<b>IT</b>	SET DI SOMMINISTRAZIONE	1	<b>HU</b>	BEADÁSI KÉSZLET	12
<b>CS</b>	SOUPRAVA PRO PODÁNÍ LÉKŮ	2	<b>LT</b>	PREPARATŲ DAVIMO RINKINYS	13
<b>DA</b>	ADMINISTRATIONSSÆT	3	<b>LV</b>	IEVADĪŠANAS KOMPLEKTS	14
<b>DE</b>	VERABREICHUNGSSET	4	<b>NL</b>	TOEDIENINGSSET	15
<b>EL</b>	ΣΕΤ ΧΟΡΗΓΗΣΗΣ	5	<b>NO</b>	ADMINISTRASJONSSETT	16
<b>EN</b>	ADMINISTRATION SET	6	<b>PL</b>	ZESTAW DO PODAWANIA	17
<b>ES</b>	SET DE ADMINISTRACIÓN	7	<b>PT</b>	CONJUNTO DE ADMINISTRAÇÃO	18
<b>FI</b>	ANTOSARJA	8	<b>RO</b>	SET DE ADMINISTRARE	19
<b>FR</b>	SETS DE PERFUSION	9	<b>SK</b>	SÚPRAVA NA PODÁVANIE LIEKOV	20
<b>HE</b>	ערכת מתן תרופות	10	<b>SR</b>	KOMPLET ZA ADMINISTRACIJU	21
<b>HR</b>	KOMPLET ZA DAVANJE LIJEKOVA	11	<b>SV</b>	ADMINISTRERINGSSET	22

REF:  
14081008

D DRUG

S SALINE / DILUENT

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**EN ADMINISTRATION SET**

**Product description:** single use administration IV sets for drugs, including chemotherapeutic solutions and hazardous drugs. The device creates a closed circuit for the safe administration of drugs.

**Intended use:** IV Sets for drugs infusion are non-invasive medical devices intended to be used for channelling of liquids for infusion and administration of drugs, including chemotherapeutic solutions, to create a customised drug solution for the needs of each patient. The IV Sets are compliant to be used in combination with Chemfort™ administration components. These devices can be connected directly or indirectly to class IIa or higher active medical devices.

**Indication for use:** The administration IV sets are designed for use with devices intended for channelling, mixing and infusing tumour treatment drugs to the patient. Chemfort™ components, included in the IV sets, prevents the introduction of microbial and airborne contaminants into the drug or fluid path for up to 7 days and 10 activations.

**Target patient population:** patients who require intravenous administration of chemotherapy drugs for therapeutic purposes. Suitable for adults, children (2-11 years) and adolescents (12-18 years).

**Users:** Medical devices designed to be used by qualified healthcare professionals.



**Contraindications:**

- Chemotherapy drugs can have various and additional side effects and contraindications. Refer to the Instruction for use of the chemotherapy drugs (e.g. antineoplastic drugs or Taxol, if photosensitive and non-photosensitive substances are present, etc.)
- These devices are not suitable for infants and neonatal.

**Warnings and precautions:**

- Do not use if the package is damaged and/or caps are missing or not in place.
- Open the packages immediately before use, always using an aseptic technique.
- The systems shall only be used by specialised personnel properly trained in the use and handling of antineoplastic chemotherapy drugs.
- This device may only be used under safe conditions in the applications and in the manner indicated in this instruction sheet.
- Use the instruments and protection systems required for handling antineoplastic chemotherapy drugs that are provided by the hospital.
- Qualified healthcare professional is responsible for choosing the modality of treatment, select the type and dosage and evaluates the risks and benefits based on the patient's condition.
- During use of the device, if any leak occurs at any point in the line or upstream of the line, or the device doesn't work properly, stop the use and replace the entire system if necessary.
- Disposable device; discard after each use. Do not reuse and/or re-sterilize. The risks caused by reuse are those of cross contamination and consequent possible infection and loss of device functionality.
- The user must follow the current legal regulations for disposal of hospital waste and hazardous drugs
- Replace tubing IV set per European Centre for Disease Prevention and Control guidelines/ local regulations and institutional protocol.

**User notice:** In the event of any serious accident related to the use of the device, this shall be reported to SIDAM and to the National Competent Authority of the Member State where the user and/or patient is established.

SYMBOLS	DESCRIPTION	SYMBOLS	DESCRIPTION
	CATALOGUE NUMBER		MEDICAL DEVICE
	LOT		CE MARKING
	DATE OF PRODUCTION		UNIQUE DEVICE IDENTIFICATION
	MANUFACTURER		QUANTITY PER BOX
	TO BE USED BEFORE		DO NOT USE IF PACKAGE IS DAMAGED
	PRODUCT STERILIZED WITH ETHYLENE OXIDE		FRAGILE, HANDLE WITH CARE
	DO NOT REUSE		KEEP DRY
	DO NOT RESTERILIZE		KEEP AWAY FROM SUNLIGHT
	LIQUID FILTER WITH PORE SIZE		KEEP IN VERTICAL POSITION
	IT DOES NOT CONTAIN NATURAL RUBBER LATEX		AUTHORISED REPRESENTATIVE IN SWITZERLAND
	IT DOES NOT CONTAIN PHTHALATES		NON-PYROGENIC
	REFER TO THE INSTRUCTIONS FOR USE		DISTRIBUTOR
	SINGLE STERILE BARRIER SYSTEM		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE

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